


Bi-Elevation[™]

PENTACAM[®] CORNEA SCLERAL PROFILE (CSP)
SOFTWARE FITTING GUIDE

Bi-Elevation[™] is an advanced scleral alignment design feature used to help reduce instances of scleral misalignment




1 VERIFY THE LENS GEOMETRY AND DIAMETER

1A. Determine lens diameter

- If HVID ≤ 11.7 mm, choose the 14.8-mm or 16.0-mm diameter lens
- If HVID > 11.7 mm, choose the 15.4-mm or 17.0-mm diameter lens

1B. Choose lens shape

- Choose **prolate** for eyes with central elevations
- Choose **oblate** for eyes with peripheral elevations or central flattening



2

DETERMINE THE AMOUNT OF BI-ELEVATION™ TO APPLY TO THE LENS

2A. Set “Ring Diameter” to the initial landing point of your desired lens diameter (in this case 13.6 mm)

14.8 = 12.2mm
15.4 = 12.8mm
16.0 = 12.8mm
17.0 = 13.6mm

2B. The “Sag. Height Astig.” value will determine the amount of Bi-Elevation™ to add to your lens SAG

3

SELECT YOUR FINAL SAG AND APS VALUES

3A. Set “Ring Diameter” to overall diameter of the lens (in this case 17.0 mm)

- Add the Bi-Elevation™ amount from step **2B** to this value for the second meridian SAG

3B. Use the remaining difference between the “Sag Height Astig.” in **3B** and the “Sag Height Astig.” in **2B** to determine the amount of toricity to add to the APS

Questions? Our expert consultants can help. Gain individualized support from our expert fitting consultants available Monday to Friday, 8:00 AM to 7:00 PM EST. Call (800) 253-3669 or email svp.consultation@Bausch.com

Visit bauschsvp.com for Important Safety Information.

®/™ are trademark of Bausch & Lomb Incorporated or its affiliates.
Any other products/brand names and/or logos are trademarks of the respective owners.
©2023 Bausch & Lomb Incorporated or its affiliates. ALZN.0034.USA.23

BAUSCH + LOMB
Zenlens
scleral lenses